
P O L I S P E C . C O M I T P H O T O N I C S . C O M



## PRODUCT CATALOGUE

REFERENCE INSTRUMENTS FOR SPECTROPHOTOMETRIC ANALYSIS

> Polispec Polispec industry Other products Installation kits



POLISPEC.COM





### LITE



POLISPEC

Getting insight

### ITPHOTONICS.COM

Founded in 2012, our company is dynamic and creative, **specialised** in spectrophotometry and applied electronics. We use our knowledge alongside theoretical and practical methods to design and build electronic systems and measuring instruments. Our core business is **applied spectrophotometry** in its various forms and for every field of application.

We created Polispec (Portable and On Line SPECtrophotometer), a line of compact, industrial-design spectrophotometers, studied and built both for portable use as well as for online installation and available for various wavelength ranges. They operate based on a luminous source interacting with the molecules and chemical bonds that characterise the matrix being analysed, thereby allowing both quantitative and qualitative measurements to be taken.

Polispec line instruments are designed to guarantee intensive use in all processes requiring immediate, precise measurements to manage variables and for self-control systems.



NIR







NIRe

PRESENTATION



SOFTWARE INNOVATION

VIS-NIT



# Polispec

# 01

LINE OF

7

Polispec LITE is an extremely compact reflection spectrophotometer, equipped with an internal lighting system and automatic references. The casing is made entirely of anodised aluminium, while a large heat sink located at the front removes the heat generated by the lighting system.

The entire system is designed for use in various settings and work conditions, thanks in part to the option of building it using different materials for the external surfaces and for the surfaces that come into contact with the products being analysed.





Dim.	21.6 x 21.3 x 8.5 cm ( <i>l x h x w</i> )
Weight	2.5 kg
Material	Anodised and/or painted aluminium
IP rating	IP68 (standard) / IP6X + IPX9K (upon request)

$\mathbf{O}$	Sensor C	MOS, 256 pix	els	
	Spectral range	e Availal	ble in va	
	Average nume	rical resoluti	on	
	Average HWHM optical resolution			
	Optical signal	collection	Direc	
	Power supply		with po angeab	
	Maximum abs	orbed power	15	
	<b>Type of measu</b> * may require exte			
	Measurement	geometry	Diffus	
	Measurement	references	Inter	
	<b>Communicatio</b> ** also available in			
	Source type	Replaceab	le halog	
$\bigcirc$	Controls and s	signals		

Hardware button dedicated to reference acquisition. Hardware button dedicated to measurement acquisition. Light and sound signal indicating acquisition status.

Software \_\_\_\_

compatible with SensoLogic and UCal Chemometric suites.

**poliPROCESS** Process analysis software compatible with SensoLogic and UCal Chemometric suites, interface available for ISObus and PLC systems (modbus over TCP/IP or S7).

01 Polispec

LITE

various ranges between 580-1080 nm

< 2 nm

< 8 nm (in range 640-1050 nm)

ct optical coupling

ower supply included,

ble rechargeable battery

5 W

ctance / transmittance \*

use / 0°

rnal and automatic

Standard RS422 nd RS485

ogen lamp

**poliDATA** Data acquisition and chemometric prediction software



**Polispec NIR is a sturdy, compact spectrophotometer with integrated reflectance measurement optics**. Designed for both manual use and in-process installation, it is built with special technological features that make it high-performing in terms of sensitivity, operating dynamics, and signal cleanliness.

These qualities make **Polispec NIR** suitable to analyse very diverse matrices, from the most reflective to the most absorbent, without requiring several versions of the same instrument for each field of application.

polispecie

Dim.	21.6 x 21.3 x 8.5 cm ( <i>l x h x w</i> )
Weight	3.2 kg
Material	Anodised and/or painted aluminium
IP rating	IP68 (standard) / IP6X + IPX9K (upon request)

Cooling control sy	vstem	Feedba
Spectrum range	900-1	700 nm
Average numerica	l resoluti	on
Average HWHM of	otical res	olution
Optical signal col	lection	Direc
Power supply		with pow angeable
 Maximum absorbe	ed power	20
<b>Type of measuren</b> * may require external a		Reflecta
Measurement geo	metry	Diffuse
Measurement refe	erences	Intern
<b>Communication cl</b> ** also available in ethe		
Source type ***	Replac	eable ha

Sensor

**Source type \*\*\*** Replaceable halogen lamp \*\*\* internal back-up lamp option available

Controls and signals

Hardware button dedicated to reference acquisition. Hardware button dedicated to measurement acquisition. Light and sound signal indicating acquisition status.

Software

**poliDATA** Data acquisition and chemometric prediction software compatible with SensoLogic and UCal Chemometric suites.

**poliPROCESS** Process analysis software compatible with SensoLogic and UCal Chemometric suites, interface available for ISObus and PLC systems (modbus over TCP/IP or S7).

10

01 Polispec

NIR

01.02

InGaAs, 256 pixels, with single-stage Peltier cooling

oack, stability T < 0.01 K

3.2 nm

3.25 nm

ect optical coupling

wer supply included, le rechargeable battery

W (

.....

tance / transmittance \*

se / 0°

nal and automatic

Standard RS422 d RS485

.....



Polispec NIRe is a sturdy, compact, extended spectral range spectrophotometer with integrated reflectance measurement optics. Designed for both manual use and in-process installation, it is built with special technological features that make it high-performing in terms of sensitivity, operating dynamics, and signal cleanliness. These qualities make **Polispec NIRe** suitable to analyse very diverse matrices, from the

most reflective to the most absorbent, without requiring several versions of the same instrument for each field of application.

01 Polispec 01.03 NIRe





Dim.	21.6 x 21.3 x 8.5 cm ( <i>l x h x w</i> )
Weight	3.3 kg
Material	Anodised and/or painted aluminium
IP rating	IP68 (standard) / IP6X + IPX9K (upon request)

\$	Sensor	-	As, 512 pixels, tage Peltier coo
	Cooling c	ontrol system	Feedback, st
	Spectrum	range 930-	2180 nm
	Average n	umerical resolut	tion 2.4 nr
	Average H	IWHM optical res	solution 4 n
	Optical si	gnal collection	Direct optica
	Power su		with power supp angeable recha
	Maximum	absorbed power	24 W
~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~	••	easurements e external accessories	
	Measuren	nent geometry	Diffuse / 0°
	Measuren	nent references	Internal and
		cation channels able in ethernet config	
		<b>pe ***</b> Repla back-up lamp option a	•
~			

 $\bigcirc$ **Controls and signals** Hardware button dedicated to reference acquisition. Hardware button dedicated to measurement acquisition. Light and sound signal indicating acquisition status.

Software \_\_\_\_

(modbus over TCP/IP or S7).

er cooling

ack, stability T < 0.03 K

2.4 nm

4 nm

optical coupling

supply included, rechargeable battery

ance / transmittance \*

al and automatic

tandard RS422 RS485

logen lamp

**poliDATA** Data acquisition and chemometric prediction software compatible with SensoLogic and UCal Chemometric suites.

poliPROCESS Process analysis software compatible with SensoLogic and UCal Chemometric suites, interface available for ISObus and PLC systems



Polispec VIS-NIT is an extremely sturdy spectrophometric sensor, designed to be able to connect to a variety of measurement probes and lighting sources via fibre optics. Being particularly compact, Polispec VIS-NIT is easy to integrate anywhere along a production line. The casing is made of anodised aluminium and is available with various levels of protection from liquids and dust. Fibre optics input with integrated optical modulator, shutter, and "neutral density" references.

CMOS, 256 pixels Sensor **A**: 340-750 nm Spectral range **B**: 580-1080 nm Average numerical resolution < 2 nm **A**: < 7 nm Average HWHM optical resolution **B**: < 10 nm **Optical signal collection** Collimated SMA input 12 Vdc with power supply included Power supply Maximum absorbed power 10 W Measurement geometry Multi-level ND attenuation. {@} Continuous optical attenuation via input beam reduction along the X axis or in predefined steps (contextual attenuation along X and Y) **Reference** acquisition Automatic internal references **Communication channel** USB

Software -

> **poliPROCESS** For in-line measurements to interface with PLC automation. Available for continuous acquisitions for filtering and for real-time chemometric predictions.

**poliDATA** Completely automated software for *stand-alone* and laboratory acquisitions and chemometric predictions.

**poliADK** Windows ADK for software developers. ADK is able to independently manage optimal instrument settings, including calculating optimal integration time. The reference acquisition procedure is also automated and is used to apply pre-treatments such as spectra smoothing or interpolation. It also extracts raw and post-processed spectral information.

Our software is compatible with SensoLogic and UCal chemometric calibration suites.

0

polrspecte

**Dim.** 110 x 10.5 x 19.5 cm (*l x h x w*)

Weight 750 g

Anodised aluminium Material

pol/spec===

14

0

: .

01 Polispec


01.04

**VIS-NIT** 



### SP3 MANAGER

### poliDATA



# 01 Polispec

### SOFTWARE

This is our main software for portable use of Polispec instruments and it offers an interface designed for use on both tablet and desktop. It is available in various languages and with two prediction engines, SensoLogic and UCal. Its basic functions are:

- Quantitative chemometric analysis: by choosing a rented and built chemometric model, a product can be analysed with different sampling methods (single swipe, multiple swipe, multipoint). The results of the analysis are instantly visible, can be exported in various formats (DAT, CSV/XML, compatible with different formulation software), and can be saved in an archive

- Analysis archive: the generated and saved data can be kept and consulted, reports can be printed, and analyses on the same product type can be compared

- Spectra acquisition: the acquired spectra can be saved to create a dataset, with export formats such as SP3 (proprietary binary software) or CSV

- Instrument diagnostics: it contains two functions for instrument diagnostics,

→ Check cell, a guided procedure that allows the user to check proper spectral calibration on the instrument

→ Diagnostic tool, a guided procedure that allows the user to make sure the main components of the instruments work (buttons, fans, internal electronics)

### СММ

### (Chemometric Model Manager)

(poliDATA). It can only be operated via a certificate that establishes the user and the accessory functions (expiry date, global standardisation...). This software has the exclusive feature of being able to use a method for "global standard**isation**", which is a special algorithm we have developed to replicate using the same calibration model on multiple instruments without having to adapt the instrument or the calibration curve itself.

This is software to import spectra saved in SP3 (or CSV) format and process them. The free version allows users to open SP3 files and convert them to CSV format to be used in conventional spreadsheets.

The licenced version has additional functions available:

- rename samples

- create an average of the samples with

the same name

- export in various formats

(CSV, CSV for Matlab or Octave, DAT, CPF)

This software also has an analysis function available (which requires a specific licence): the selected spectra can be predicted using a chemometric model, the analyses can be exported in CSV format or copied/pasted directly into a spreadsheet. There is the option of whether or not to display the statistical parameters of the analyses for each constituent.

Finally, there is a **dataset function** available (which requires a specific licence): by opening an SP3 file and adding the analyses (copying from a spreadsheet or importing the file in CSV format), a dataset file can be created and exported in DAT or CPF (SensoLogic) format.

This software allows users to prepare chemometric models that can be used by the previously described software


# Polispec industry

# 02

LINE OF

Polispec LITE industry is an extremely compact reflection spectrophotometer, equipped with an internal lighting system and automatic references. Especially suitable for macroelement analysis, **Polispec LITE industry** easily adapts to various working conditions. The casing is made entirely of anodised aluminium, while a large heat sink located at the front removes the heat generated by the lighting system.





Dim.	21.6 x 21.3 x 8.5 cm ( <i>l x h x w</i> )
Weight	2.5 kg
Material	Anodised and/or painted aluminium
IP rating	IP68 (standard) / IP6X + IPX9K (upon request)

$\mathbf{O}$	Sensor	CMOS, 256 pixels		
	Spectral rai	ige	Availat	ole in va
	Average numerical resolution			
	Average HWHM optical resolution			
	Optical sign	al colle	ction	Direct
	Power supp	ly	12 Vdc v interch	with pov angeabl
	Maximum a	bsorbed	power	15
	<b>Type of mea</b> * may require e			Reflect
	Measureme	nt geom	etry	Diffus
	Measureme	nt refer	ences	Interi
	<b>Communica</b> ** also availabl			
	Source type	Re	placeab	le halog
<u>-</u>				

 $[\bigcirc$ Controls and signals

> Available in two versions, with or without processing unit integrated into the instrument case.

\_\_\_\_ Software

> **poliPROCESS** Process analysis software compatible with SensoLogic and UCal Chemometric suites, interface available for ISObus and PLC systems (modbus over TCP/IP or S7).

LITE industry arious ranges between 580-1080 nm

< 2 nm

< 8 nm (in range 640-1050 nm)

ct optical coupling

ower supply included,

ble rechargeable battery

5 W

tance / transmittance \*


use / 0°

rnal and automatic

Standard RS422 nd RS485

gen lamp





Feedback, stability T < 0.01 K

3.2 nm

Direct optical coupling


Reflectance / transmittance \*

Internal and automatic



**Polispec NIRe industry is a sturdy, compact, extended spectral range spectrophotometer with integrated reflectance measurement optics**. Designed for in-process installation, it is built with special technological features that make it high-performing in terms of sensitivity, operating dynamics, and signal cleanliness.





Dim.	21.6 x 21.3 x 8.5 cm ( <i>l x h x w</i> )
Weight	3.3 kg
Material	Anodised and/or painted aluminium
IP rating	IP68 (standard) / IP6X + IPX9K (upon request)

	uge i citic
system	Feedbac
930-:	2180 nm
al resoluti	ion 2
optical res	olution
llection	Direct o
	vith power angeable r
oed power	24 V
	Reflectar
ometry	Diffuse
ferences	Interna
	system 930-: al resolution optical res llection 12 Vdc v

Sensor

Controls and signals

Available in two versions, with or without processing unit integrated into the instrument case.

**Software** 

**poliPROCESS** Process analysis software compatible with SensoLogic and UCal Chemometric suites, interface available for ISObus and PLC systems (modbus over TCP/IP or S7).

Dual chip InGaAs, 512 pixels, with double-stage Peltier cooling ntrol system Feedback, stability T < 0.03 K range 930-2180 nm umerical resolution 2.4 nm WHM optical resolution 4 nm mal collection Direct optical coupling ply 12 Vdc with power supply included, interchangeable rechargeable battery absorbed power 24 W

e / 0°


nal and automatic

tandard RS422 <sub>RS485</sub>

alogen lamp

.....





< 2 nm

**A**: < 7 nm **B**: < 10 nm

Collimated SMA input

Multi-level ND attenuation. Continuous optical attenuation via input beam reduction along the X axis or

Automatic internal references

RS422 / RS485 / RS232 (to be determined while ordering)

**poliPROCESS** For in-line measurements to interface with PLC automation. Available for continuous acquisitions for filtering and for real-time chemom-

automated and is used to apply pre-treatments such as spectra smoothing




Maintaining the same industrial and compact design as the standard version, Polispec LITE Industry CONTACTLESS is a very versatile system for application on all processes where non-contact measurement is necessary. Adaptable to a wide range of distances from the product, this product line is equipped with its own illuminant with internal feedback for automatic references during the process. The large framed area and the optical geometry that characterize it, allow optimal work even if the relative distance between the product and the measurement optics undergoes variations.

pol/speces

Dim. 40,5 x 8,5 x 25 cm (*l x h x w*) Weight 5 kg

Anodised and/or painted aluminium Material Instrument IP68 (standard) / IP6X + IPX9K IP rating (upon request); CONTACTLESS module IP55



Software =

> poliPROCESS Process analysis software compatible with SensoLogic and UCal Chemometric suites, interface available for ISObus and PLC systems (modbus over TCP/IP or S7).

02 Polispec industry

02.05

LITE

industry

CONTACTLESS

Available in various ranges between 580-1080 nm

< 2 nm

< 8 nm (in range 640-1050 nm)

Direct optical coupling

**Power supply** Instrument: 12 Vdc with power supply included, interchangeable internal battery CONTACTLESS module: 12 Vdc 20 W with dedicated power

Instrument 20 W, CONTACTLESS module 20 W

Angular / Diffuse

Internal and automatic

Standard RS422

< 1



By exploiting the entire spectral range of Polispec NIR and borrowing the same industrial design and compactness of the standard version, **Polispec NIR Industry CONTACTLESS** is a very versatile system for application on all processes where non-contact measurement is necessary. Adaptable to a wide range of distances from the product, this product line is equipped with its own illuminant with internal feedback for automatic references during the process. The large framed area and the optical geometry that characterize it allow optimal work even if the relative distance between the product and the measurement optics undergoes variations.

# 02 Polispec industry

02.06

NIR industry

### CONTACTLESS





Dim.	40,5 x 8,5 x 25 cm ( <i>l x h x w</i> )
Weight	4,3 kg

Material Anodised and/or painted aluminium IP rating Instrument IP68 (standard) / IP6X + IPX9K (upon request); CONTACTLESS module IP55

<b>Sensor</b> InGaAs,	256 pixe	els, with si
Spectrum range	900-170	)0 nm
Average numerical r	esolutio	<b>n</b> 3.2
Average HWHM optic	cal resol	ution
Optical signal collec	tion	Direct op
Maximum absorbed	power	Instrur
Measurement geome Measurement refere Communication chan ** also available in etherne Source type Rep Working distance	etry nces nnels ** et configura placeable From	Angular / Internal Stan tion or RS4 e halogen 20 cm to 0
	Spectrum range Average numerical r Average HWHM optic Optical signal collect Power supply In able internal battery supply Maximum absorbed Masurement geome Measurement geome Measurement refere Communication chan ** also available in etherne Source type Rep	Maximum absorbed power Type of measurements Measurement geometry Measurement references Communication channels ** ** also available in ethernet configura

Controls and signals  $\Theta$ 

the instrument case

Software =

(modbus over TCP/IP or S7).

single-stage Peltier cooling

2 nm

3.25 nm

optical coupling

with power supply included, interchangenodule: 12 Vdc 20 W with dedicated power

ument 20 W, CONTACTLESS module 20 W

nce

/ Diffuse

l and automatic

ndard RS422 6485

lamp

60 cm

< 1

Available in two versions, with or without processing unit integrated into

**poliPROCESS** Process analysis software compatible with SensoLogic and UCal Chemometric suites, interface available for ISObus and PLC systems



### poliPROCESS

### SP3 MANAGER

02 Polispec industry

### SOFTWARE

### This is the main software to use instruments on processes

(industrial plants or machinery) requiring continuous product analysis. The software is modular and can be configured or expanded (by custom-developing new modules) to adapt to specific needs. The main reference modules are:

### GPS module

which acquires NMEA data from GPS antennas via serial connection and uses the position received to geo-reference the analyses

### Datalogger module

which saves the analyses in a file in CSV or KML (Keyhole Markup Language) formats

### ISObus module

which manages the analyses via an application on the Virtual Terminal (VT); the acquired data are saved from the device task controller and can be exported in ISOXML file and used in most analysis platforms. The VT interface can be used to control the process completely: selecting the product to analyse (chemometric model to use), displaying instantaneous and average analysis (of the process), receiving notifications and alarms for any problems

### PLC module

which controls the measurement process via PLC using two different protocols: modbus over TCP/IP or S7 and ISO over TCP/IP or S7 (Siemens S7). The measurement process reports the analysis values and the system status on registers that can be configured during installation


### СММ

(Chemometric Model Manager)









This is software to import spectra saved in SP3 (or CSV) format and process them. The free version allows users to open SP3 files and convert them to CSV format to be used in conventional spreadsheets.

The licenced version has additional functions available:

- rename samples

- create an average of the samples with

the same name

- export in various formats

(CSV, CSV for Matlab or Octave, DAT, CPF)

This software also has an analysis function available (which requires a specific licence): the selected spectra can be predicted using a chemometric model, the analyses can be exported in CSV format or copied/pasted directly into a spreadsheet. There is the option of whether or not to display the statistical parameters of the analyses for each constituent.

Finally, there is a **dataset function** available (which requires a specific licence): by opening an SP3 file and adding the analyses (copying from a spreadsheet or importing the file in CSV format), a dataset file can be created and exported in DAT or CPF (SensoLogic) format.

This software allows users to prepare chemometric models that can be used by the previously described software (poliPROCESS). It can only be operated via a certificate that

establishes the user and the accessory functions (expiry date, global standardisation...).

This software has the exclusive feature of being able to use a method for "global standardisation", which is a special algorithm we have developed to replicate use of the same calibration model on multiple instruments without having to adapt the instrument or the calibration curve itself.





Focus

### PROBE

for: Polispec Polispec industry All the instruments in the Polispec and Polispec industry lines can be fitted with different scanning window configurations in order to make the system compatible for manual use, in-line installation, or mixed use. Depending on what is needed, the instruments can be configured with the following types of scanning windows: - flat (0 mm) - extended probe (2.5 mm, 11 mm, 14 mm, 17 mm and custom) The probes, with a standard diameter of 50 mm, are usually


fitted with a sapphire scanning window and can be made from various materials, depending on the application requirements.

MECHANICALLY ATTACHED S C A N N I N G WINDOW





In the standard configuration, the sapphire scanning window is glued directly onto the probe with special glues that ensure resistance and adherence.



However, for more difficult settings or where, for legislative reasons, a glued-on window cannot be used, there is a special probe available where the scanning window is fastened mechanically.

Depending on the requirements, the probes can be made from various materials (i.e., food grade POM, stainless steel, hardox) in order to ensure the system is compatible in terms of food safety, thermal insulation, and resistance to corrosion or wear.

# Polispec

# 03

O T H E R P R O D U C T S , A C C E S S O R I E S A N D E X T E N S I O N K I T S

## HARDWARE AND SOFTWARE SYSTEMS TO INTERFACE POLISPEC INSTRUMENTS WITH OPERATING MACHINERY

### 03.03





**POWER BOX** is a power supply control unit that ensures the Polispec instruments function properly when they are installed on self-propelled or towing machinery. POWER BOX manages proper current flow during system operation, as well as allowing the sensor to switch on and off via the key or engine on/off signal.

# 03

OTHER PRODUCTS, ACCESSORIES AND EXTENSION KITS



Transferring data from the ISObus Task Controller to the cloud is no longer a problem, thanks to UBeQuo.

Simply connect this small device to the machine terminal USB port and connect it to WiFi or to your mobile hotspot network to have the data immediately available on the cloud and visible directly from the office, regardless of where the machine is located.



AGS (Active Grain Sampler) is an accessory designed to install Polispec systems on combines or grain elevators. It operates to ensure continuous, representative product flow for the corresponding measurement, which self-regulates based on the instantaneous flow of material in the main duct.

### 03.02





1 UBeQUE

ISO BOX is an intelligent control unit that, in addition to managing the power supply to the Polispec sensors when installed on self-propelled or towing machinery, processes the signals it receives from the sensor and integrates sensor operation in the machine ISObus network. The same control unit can be used to integrate systems in proprietary CANbus networks.

OTHER PRODUCTS

M U L T I - P U R P O S E " A D D - O N " A C C E S S O R I E S



03.18

## DESKTOP TABLE WITH ROTATIONG SAMPLING CELL

This support transforms Polispec instruments into desktop measurement systems. It is a simple structure with a modern design that hold the spectrophotometer vertically to be used as a real benchtop system. A practical measuring cell allows to contain the sample while it is rotating for the measure. Once the measurement is completed, the cell is removed allowing easy cleaning of the plate also thanks to the special slide discharge. 03.06

### PROTECTIVE PLATE KIT

The base of this magnetically attachable accessory is wider than the instrument itself and extends with two practical folds on the front and back of the instrument in order to protect it from stubborn dirt during use. The accessory is designed to allow the system to be cleaned quickly when it is used on muddy or slimy matrices.



03.19

# DESKTOP BASE FOR VERTICAL POSITIONING

This desktop accessory makes it easy to quickly position the instrument vertically to measure small samples, cells, Petri dishes, or to use other accessories.

" A D D - O N "
FILLER ELEMENT
This filler element easily attaches to the bottom of the instrument magnetically and allows instruments with an extended probe installed on them to be used portably.
The accessory is compatible with the Polispec LITE, Polispec NIR and Polispec NIRe

The accessory is compatible with the **Po** instruments.

03.07





03.16

## CONVEX SAMPLE CONTACT PROBE WITH REDUCTION OF THE FRAMED AREA

This magnetically attachable accessory allows the instrument's framed area to be reduced to an area of Ø5 mm or Ø10 mm. Very practical to use, it also lends itself well to measuring convex surfaces like, for example, small fruits and solid objects with non-planar surfaces.

## RETURN PROBE FOR 90° POSITIONING

This element returns both the light beam irradiating the sample and the collected signal at 90°, allowing the instrument to be positioned horizontally compared to its original optical geometry.

THE PRACTICAL SCREW SYSTEM ALLOWS THE LENS TO BE REPLACED QUICKLY TO CHANGE THE FRAMED **AREA** (Ø5 mm, Ø10 mm AND CUSTOM Ø AVAILABLE)









03.17



### MULTI-PURPOSE "ADD-ON" ACCESSORIES



# "ADD-ON" KIT FOR LIQUID READING

This accessory, which is compatible with **Polispec LITE**, **Polispec NIR** and **Polispec** NIRe instruments in portable configuration, consists of a magnetic cuvette-holder support that attaches to the front of the instrument's scanning window and measures liquids, both in reflectance and transflectance, as it fits a reflector that can be removed from the cuvette and positioned at a variable optical path.

THE CUVETTE, ADJUSTABLE ON BOTH SIDES, ALLOWS THE REFLECTOR TO BE POSITIONED AT A VARIABLE **OPTICAL PATH**, WHICH CAN BE INCREASED WITH 2 mm PITCH



## THE CUPIS 80 mm TALL AND IS **SUITABLE BOTH** FOR INSTRUMENTS WITH HØ PROBES AND INSTRUMENTS WITH EXTENDED PROBES



03.20

# CUP FOR LIQUID READING WITH DRIP TRAY

This accessory is used to analyse liquid samples of varying nature. A holding cup, located inside a tray, contains liquid that makes use of the "meniscus effect" to adhere to the scanning window perfectly. Any excess liquid that could spill out is collected by the surrounding tray that can be easily removed, emptied, and washed.







03.21

# KIT FOR MEASUREMENT OF LIQUIDS IN TEST TUBE

This practical accessory allows the measurement of small quantities of liquid inserted into a washable glass test tube. The Polispec device is placed vertically so as to accommodate a magnetic support that holds the test tube. The measurement is then carried out in transflectance, inserting a special golden reflector immersed in the liquid. The optical path of the reflector can be adapted to specific needs.

# P O R T A B L E K I T



# WIRING KIT

03.13

for installation on process/machinery

# TABLET VEHICLE DOCKING STATION KIT

03.14

tablet + vehicle docking station + high density cable

# P R O C E S S K I T

### 03.12

AC 100/240 V to DC 12 V power supply + case + WiFi antenna (in the event of standard instrument) + USB cable + poliPROCESS SW



# TABLET VESA DOCKING STATION KIT

03.15

tablet + VESA docking station







49

# Polispec

# 04

INSTALLATION KIT



# 04

# INSTALLATION KIT for:


Polispec Lite Polispec NIR Polispec NIRe All instruments in the Polispec line have various installation kits available to guarantee proper system installation, based on where they are situated, on operating machinery or on industrial production lines.

Their purpose is to guarantee safe, proper installation and, depending on the requirements, to make it easy to remove the sensor for cleaning, diagnostics, or use as a portable system. This kit consists of an **mounting system with brackets to install Polispec instruments in-process or on operating machinery** (for example, a forage harvester) and in situations requiring horizontal installations that involve quickly removing the instrument for inspection or for use as a portable system.



INSTALLATION KIT





53







Mounting system complete with adapter flange to weld, available in various thicknesses (6 mm standard). Recommended when installing on surfaces that are not flat or insufficient to sustain direct installation and in cases requiring the instrument to potentially be removed for portable use.




Kit for on-duct installation, available in various formats and various flanging (6" and 8" standard). The kit is supplied complete with sealing gaskets and the accessories needed to install and protect the **Polispec systems**. The standard kit is suitable for operating pressure up to 3 BAR, with pressure peaks up to 7.5 BAR.

INSTALLATION KIT



ΚΙΤ

### OTHER INSTALLATIONS















# OPTIONAL



Protective cover with optional forced ventilation, designed for settings that must ensure forced air circulation to better dissipate the heat generated by the instrument  $^{\circ}T > 40-45^{\circ}C$ 






INSTALLATION KIT 04.06

ΚΙΤ

This flat flange universal coupling system, compatible with instruments in the Polispec Lite, Polispec NIR and Polispec NIRe lines, makes it possible to install the systems even in the tightest spaces. This flat flange can be adjusted to every position and is attached at the base of the instrument, which allows the system to be anchored using the four through-holes for M5 screws and short "wings" that protrude sideways with reference to the instrument.

THE SYSTEM MAKES IT EASY TO INSTALL ON SQUARE PIPING, HOPPERS, OR FLAT SURFACES REGARDLESS OF THE ORIENTATION (VERTICAL, HORIZONTAL, LATERAL)



This installation kit for 90° positioning, suitable for **Polispec Lite**, **Polispec** NIR and Polispec NIRe instruments, allows the system to be positioned horizontally and requires a return probe for 90° positioning (03.17). The kit is equipped with a protective box in white plastic and is designed to allow the system to be removed quickly for cleaning and maintenance.

KIT



ΚΙΤ

# S E C T O R S O F A P P L I C A T I O N

Our passion and our specialized expertise allow us to create **targeted and effective solutions** to meet the needs, even the most particular ones, in the field of spectrophotometric analysis, technological integration and electronic and software applications. We are able to reach even the most ambitious objectives, regardless of the application sector, thanks to the distinct professionalism that distinguishes us and the recognized quality of our technology.



INDUSTRIAL



A G R I C U L T U R A L A N D A G R I - F O O D



### CHEMICAL AND SCIENTIFIC



POLISPEC.COM

ITPhotonics S.r.l. via Astico 39, 36030 Fara Vicentino (VI) - Italy T. (+39) 0445-1925221 M. info@itphotonics.com

### I T P H O T O N I C S . C O M